



# MASTIDIAN JERINGAS

## INTRAMAMMARY SUSPENSION

One of the infectious processes which have a major economic and sanitary impact in dairy cattle are mammary gland infections and, more specifically, inflammatory infections known as **MASTITIS**.

Mastitis is extremely common and refers to an inflammation of the udder as a response to infection. These infections are caused mainly by Gram (+) and Gram (-) microorganisms, each group acting separately or together.

Lactating dairy cows and dry cows are both susceptible to mastitis, but the disease is more common in lactating cows because of greater susceptibility and because it is transmitted primarily at milking. Mastitis is also more apparent in lactating cows due to its impact on milk production.

Once a process of mastitis occurs, there is no alternative other than to administer a treatment which directly controls the infectious agents.

### COMPOSITION

Per 1 syringe.

- ◆ Cloxacillin sodium .....200 mg
- ◆ Colistin (sulphate) ..... 400.000 I.U.
- ◆ Excipient to ..... 10 ml

### DOSAGE AND ADMINISTRATION ROUTE

Intramammary route.

After milking thoroughly, apply the contents of one syringe in each affected quarter. At the first administration is advisable to treat all quarters, even if some of them are not affected. Repeat the administration after 24 hours, for 3-5 days.

### PROPERTIES

#### MASTIDIAN JERINGAS

consists of a synergical association of two bactericidal antibiotics: Cloxacillin, active against Gram(+) microorganisms, and Colistin which acts against Gram(-).

The result of this association is an effective product which offers broad-spectrum therapy for lactating cows being approved to treat CATTLE mastitis caused by *Staphylococcus*, *Streptococcus*, *Enterobacteria*, *Pseudomonas* and *Corynebacterium*.



### WITHDRAWAL PERIOD

MEAT: 4 days  
MILK: 3 days

### PRESENTATION

- . Boxes with 4 and 12 single-dose syringes.
- . Drums with 100 single-dose syringes.



**SUPER'S DIANA, S.L.**

Ctra. C-17, km 17  
08150 Parets del Vallès (BARCELONA) SPAIN  
Tel. +34 93 8437536 / Fax. +34 93 8436545  
E-mail: [export@supersdiana.com](mailto:export@supersdiana.com)  
[www.supersdiana.com](http://www.supersdiana.com)